

# Therapy Comparison

Key

Feature ■ Safe ■ Caution ■

Feature ■ Safe ■ Caution ■

Feature ■ Safe ■ Caution ■

	PEMF		
	Feature	Safe	Caution
Pain Relief Mode	<span style="color: blue;">■</span>		
Wound Relief Mode	<span style="color: blue;">■</span>		
Stress Relief Mode	<span style="color: blue;">■</span>		
Need to use Conductor Gel			
Gel Allergy potential			
Fully portable	<span style="color: blue;">■</span>		
Use on Acupuncture points	<span style="color: blue;">■</span>	<span style="color: green;">■</span>	
Use with Prosthetic Implants	<span style="color: blue;">■</span>	<span style="color: green;">■</span>	
Contraindicate: Chemo Treatment			<span style="color: red;">■</span>
Contraindicate: Menieres Disease		<span style="color: green;">■</span>	
Usable on Skull area [1]	<span style="color: blue;">■</span>	<span style="color: green;">■</span>	
Potential bone scorch [2]		<span style="color: green;">■</span>	
Usable on Reproductive Organs [3]	<span style="color: blue;">■</span>	<span style="color: green;">■</span>	
Usable for treating Wounds [4]	<span style="color: blue;">■</span>	<span style="color: green;">■</span>	
Usable with cardiac pacemaker [5]			<span style="color: red;">■</span>
Usable on Varicose Veins [6]	<span style="color: blue;">■</span>	<span style="color: green;">■</span>	
Usable Epiphysis growth plates [7]	<span style="color: blue;">■</span>	<span style="color: green;">■</span>	
Clinical Studies on Therapy	<span style="color: blue;">■</span>		
Clinical Study on Devices	<span style="color: blue;">■</span>		

	Ultrasound		
	Feature	Safe	Danger
Pain Relief Mode	<span style="color: blue;">■</span>		
Wound Relief Mode			
Stress Relief Mode			
Need to use Conductor Gel	<span style="color: blue;">■</span>		
Gel Allergy potential			<span style="color: red;">■</span>
Fully portable			
Use on Acupuncture points			<span style="color: red;">■</span>
Use with Prosthetic Implants			<span style="color: red;">■</span>
Contraindicate: Chemo Treatment			<span style="color: red;">■</span>
Contraindicate: Menieres Disease			<span style="color: red;">■</span>
Usable on Skull area [1]			<span style="color: red;">■</span>
Potential bone scorch [2]			<span style="color: red;">■</span>
Usable on Reproductive Organs [3]			<span style="color: red;">■</span>
Usable for treating Wounds [4]			<span style="color: red;">■</span>
Usable with cardiac pacemaker [5]			<span style="color: red;">■</span>
Usable on Varicose Veins [6]			<span style="color: red;">■</span>
Usable Epiphysis growth plates [7]			<span style="color: red;">■</span>
Clinical Studies on Therapy	<span style="color: blue;">■</span>		
Clinical Study on Devices			

	TENS		
	Feature	Safe	Danger
Pain Relief Mode	<span style="color: blue;">■</span>		
Wound Relief Mode			
Stress Relief Mode			
Need to use Conductor Gel	<span style="color: blue;">■</span>		
Gel Allergy potential			<span style="color: red;">■</span>
Fully portable			
Use on Acupuncture points			
Use with Prosthetic Implants			<span style="color: red;">■</span>
Contraindicate: Chemo Treatment			<span style="color: red;">■</span>
Contraindicate: Menieres Disease			<span style="color: red;">■</span>
Usable on Skull area [1]			<span style="color: red;">■</span>
Potential bone scorch [2]			
Usable on Reproductive Organs [3]			<span style="color: red;">■</span>
Usable for treating Wounds [4]			
Usable with cardiac pacemaker [5]			<span style="color: red;">■</span>
Usable on Varicose Veins [6]			
Usable Epiphysis growth plates [7]			
Clinical Studies on Therapy	<span style="color: blue;">■</span>		
Clinical Study on Devices	<span style="color: blue;">■</span>		

[1] Ultrasound should not be applied to the brain, spinal cord or large subcutaneous peripheral nerves.

[2] If any pain or uncomfortable "prickly" sensation is felt by the patient, this may be an indication that the bones or nerve endings in the vicinity of the ultrasonic beam are becoming, or are already, overheated. In this case the ultrasound power or intensity must be immediately reduced.

[3] Ultrasound physical therapy should not be applied to the reproductive organs.

[4] Treatment of acute infection of bone or soft tissue should not be carried out as the treatment could force areas of pus into surrounding tissue, thereby spreading infection.

[5] Ultrasound should not be used in the thoracic area if the patient has a cardiac pacemaker of any kind. PEMF can be used 20cm or more away from site.

[6] Blood vessels in poor condition should not be treated with Ultrasound as the vessel walls may rupture as a result of the exposure.

[7] Ultrasound should not be used directly over the active Epiphysis regions (growth plates), which are in every major joint. This is even more important where children are concerned.